

Stabilab™ - the new Standard in Tartrate Stability testing

Rodney Marsh
Managing Director
Winechek
Busselton
Western Australia
rodney@winechek.com

Introduction

In recent years many advances in tartrate stabilization have occurred giving winemakers more options than the expensive cold contact process. These include electrodialysis, the use of mannoproteins (Mannostab™, Laffort Australia) and shortly CMC carboxymethylcellulose which is currently with ANZFA for approval in Australia. It would seem that there might also be need to change and upgrade our testing procedures from the old refrigeration/freezer test.

At Winechek we have found this to be true and in February 2010 we installed the latest technology, a Stabilab™ (Figure 1) machine from France.

Figure 1. The Stabilab™ with conductivity probes & agitator assembly in place



Background

At Winechek we gave up on the traditional and widely used refrigeration or freezer tests in 2005 as we found them neither reliable, nor easily repeatable. Our findings are confirmed by P.Ribereau-Gayon et al (2006) & B.W. Zoeckli et al (1997) and many others in the documented literature as the test is based on the phenomenon of spontaneous non-induced crystallization, a slow undependable/unpredictable process.

We have mainly used the Saturation Temperature (Tsat) conductivity test established by Erbsloh to estimate Stabilization Temperature (Tcs) for tartrate stability testing. The Tsat test has the advantage of being rapid, 5-10 minutes test time and is therefore very useful when a winery is utilizing the Cold contact process to tartrate stabilize their wine. We found that it was excellent for tank whites, which is the majority of samples but possible gave false failing results with some dark Roses, most reds and extended lees contact whites. We investigated other means of determining tartrate stability including Potassium measurement pre & post Potassium Hydrogen Tartrate (KHT) seeding and chilling, also the -4°C for 3 & 6 days but found all these methods either costly, hard to interpret or too long to be useful for our clients. We determined to find another rapid method that would fill this whole in our testing ability and after about a year of investigation finally found the Stabilab™.

The Stabilab™ was researched and developed by INRA (one of the largest agricultural research institutes in Europe) & Eurodia (producers of the Electrodialysis technology) to model tartrate stability of wine at -4°C for 6 days, the OIV/European reference test. It is basically a strictly temperature controlled cold bath with agitators to constantly stir the wine under test & conductivity probes (Figure 2) to constantly measure changes in electrical conductivity ($\mu\text{S}/\text{cm}$). The conductivity of a wine is an objective measure of the concentration of ions in solution, Potassium (K) is the major cation found in wine normally more than ten fold the ion concentrations of Sodium & Calcium. When the precipitation of Potassium Hydrogen Tartrate (KHT) crystals occurs we observe a decrease in Potassium & therefore conductivity. By supersaturating the wine with KHT crystals, holding the wine at -4°C and continuous agitating we optimize & facilitate the precipitation of tartrate crystals.

Figure 2. The conductivity probe & agitator assembly, 8 sets for 8 tests.



The Stabilab™ Tests

The Stabilab™ can perform two different tests, **DIT%** ("Degré d'Instabilité Tartrique" or Degree of Tartrate Instability test) and the **ISTC 50** (*Indice de Stabilité Tartarique Critique*, or Critical Tartaric Stability Index with a 50g/hL KHT addition).

DIT%

The purpose of the test is to answer the following questions:

- Is the wine stable? (Does not require treatment)
- If the wine is unstable, what level of treatment is required?

| | Procedure | What is happening |
|----|--|--|
| 1. | Temperature of the wine sample is measured and must be between 15° & 25°C | Getting a conductivity reading at room temperature for comparison to -4°C |
| 2. | Wine sample brought to -4°C & then another conductivity measurement | Preparing sample for actual determination of stability & getting an initial reading |
| 3. | Add a quantity equivalent to 4g/L of Potassium Hydrogen Tartrate (KHT) to the wine | Supersaturating the wine with KHT & have a large amount of seeding crystals available to induce precipitation of KHT |
| 4. | The wine is agitated at -4°C for 4 hrs with | Create a perfect environment for |

| | | |
|----|--|---|
| | conductivity measurements made every 10 minutes | the formation of KHT crystals & their precipitation. The monitoring of conductivity indicates changes in the concentration of Potassium (K) caused by the precipitation of KHT. |
| 5. | The resulting conductivities are transferred to a computer software package that calculates the DIT% and gives many other measures | Assessment of wine stability |

The software generates a report showing the results (Figure 3)

Figure 3 Example of a DIT% report

Test nouveau logiciel the 01/07/2008 at 4:26 pm



DIT


Echantillon 1

Wine presumed UNSTABLE

01/02/2008

| °C | min. | C meas. | C conv. | C calc. |
|------|------|---------|---------|---------|
| -4.0 | 0 | 918 | 1999 | 1970 |
| -4.0 | 10 | 828 | 1803 | 1783 |
| -4.0 | 20 | 783 | 1705 | 1708 |
| -4.0 | 30 | 760 | 1655 | 1667 |
| -4.0 | 40 | 750 | 1633 | 1641 |
| -4.0 | 50 | 743 | 1618 | 1624 |
| -4.0 | 60 | 738 | 1607 | 1611 |
| -4.0 | 70 | 734 | 1599 | 1601 |
| -4.0 | 80 | 732 | 1594 | 1593 |
| -4.0 | 90 | 729 | 1588 | 1587 |
| -4.0 | 100 | 727 | 1583 | 1582 |
| -4.0 | 110 | 725 | 1579 | 1578 |
| -4.0 | 120 | 723 | 1575 | 1574 |
| -4.0 | 130 | 722 | 1572 | 1571 |
| -4.0 | 140 | 721 | 1570 | 1569 |
| -4.0 | 150 | 720 | 1568 | 1566 |
| -4.0 | 160 | 719 | 1566 | 1564 |
| -4.0 | 170 | 717 | 1562 | 1562 |
| -4.0 | 180 | 716 | 1559 | 1560 |
| -4.0 | 190 | 716 | 1559 | 1559 |
| -4.0 | 200 | 716 | 1559 | 1558 |
| -4.0 | 210 | 715 | 1557 | 1556 |
| -4.0 | 220 | 714 | 1555 | 1555 |
| -4.0 | 230 | 713 | 1553 | 1554 |
| -4.0 | 240 | 712 | 1551 | 1553 |

| | |
|--------------------------------|--------|
| Initial Conductivity (µS/cm) | 1999 |
| Coefficient of conversion | 2.178 |
| Coefficient of correlation | 0.9984 |
| Conductivity after 4h (µS/cm) | 1553 |
| Conductivity drop after 4h | 22.3% |
| Conductivity after 24h (µS/cm) | 1534 |
| Conductivity drop after 24h | 23.3% |
| Infinite Conductivity (µS/cm) | 1530 |
| Rate of electro dialyse | 23.5% |
| Initial temperature (°C) | 20.0 |



This report shows us on the left side a table with 5 columns:
 - the 1st column presents the temperature in the wine during the test (the set point),

- the 2nd column presents the time of the test,
- the 3rd column presents the measured conductivity at the temperature of the test,
- the 4th column presents the converted conductivity by the coefficient of conversion (conductivity converted to 20°C),
- the 5th column presents the calculated conductivity from the mathematical model.

On the right side, calculated results are presented:

- The initial measured conductivity at room temperature (**C₀**). This conductivity is compensated at 20°C
- The coefficient of conversion Conductivity_{20°C}/Conductivity_{-4°C}.
- The coefficient of correlation between the converted curve and the calculated curve (r²),
- The conductivities from the calculated curve after 4, 24 hours and infinite time (**C_∞**),
- The DIT (“Degré d’Instabilité Tartrique” = conductivity drop) after 4, 24 hours and infinite time. The DIT% is the Rate of Electrodialysis.
- The initial temperature of the wine

The graph is a representation of the change in conductivity vs time of the test

ISTC 50

The purpose of the test is to answer the following questions:

- Is the wine stable after some stabilization treatment

| | Procedure | What is happening |
|----|---|--|
| 1. | To a quantity of wine the equivalent of 0.5g/L (50g/hL) KHT is added | Slightly destabilize the wine |
| 2. | Place wine in a 36°C waterbath for 30 minutes | Completely dissolve the KHT |
| 3. | Measure the conductivity at this temperature | Getting a conductivity reading at 36°C temperature for comparison to -4°C |
| 4. | Wine sample brought to -4°C & then another conductivity measurement | Preparing sample for actual determination of stability & getting an initial reading |
| 5. | Glass beads added | To induce precipitation of KHT by vibration |
| 6. | The wine is agitated at -4°C for 2 hrs for whites & Rose & 4hrs for Red wines with conductivity measurements made every 10 minutes. | Create a perfect environment for the formation of KHT crystals & their precipitation. The monitoring of conductivity indicates changes in the concentration of Potassium (K) caused by the precipitation of KHT. |
| 7. | The resulting conductivities are transferred to a computer software package that calculates the ISTC50 result | Assessment of wine stability |


The software generates a report showing the ISTC50 results (Figure 4)

Figure 4 Example of a ISTC50 report

| °C | min. | C meas. |
|------|------|---------|
| -4.0 | 0 | 1001 |
| -4.0 | 10 | 998 |
| -4.0 | 20 | 988 |
| -4.0 | 30 | 985 |
| -4.0 | 40 | 984 |
| -4.0 | 50 | 983 |
| -4.0 | 60 | 982 |
| -4.0 | 70 | 982 |
| -4.0 | 80 | 981 |
| -4.0 | 90 | 981 |
| -4.0 | 100 | 981 |
| -4.0 | 110 | 981 |
| -4.0 | 120 | 980 |
| -4.0 | 130 | 980 |
| -4.0 | 140 | 980 |
| -4.0 | 150 | 979 |
| -4.0 | 160 | 979 |
| -4.0 | 170 | 977 |
| -4.0 | 180 | 974 |
| -4.0 | 190 | 971 |
| -4.0 | 200 | 967 |
| -4.0 | 210 | 958 |
| -4.0 | 220 | 946 |
| -4.0 | 230 | 927 |
| -4.0 | 240 | 901 |

Kind of wine : Red

Wine presumed UNSTABLE



In the left hand table;
 First column is temperature in °C
 Second column is time in minutes
 Third column is conductivity readings in uS/cm

The graph is a representation of the change in conductivity vs the time of the test.


Results

The wine is STABLE if the conductivity drop is $< 3 \mu\text{S}/\text{cm}$ between 40 and 120 minutes for White and Rosé wines.

The wine is STABLE if the conductivity drop is $< 5 \mu\text{S}/\text{cm}$ between 40 minutes and 240 minutes for Red wine.

Obviously the wine in Figure 4 was found unstable with a conductivity drop of $83 \mu\text{S}/\text{cm}$ between 40 and 240 minutes. Stable wine graphs appear like Figure 5. Note the lack of conductivity change for most of the time of testing.

Figure 5 ISTC50 graph of stable wine



Interestingly at the end of the ISTC50 test when cleaning up the machine and disposing of the sample it is very comforting to observe the tartrate precipitates or lack thereof which correlates with the results see Image 1. The sample on the left is very unstable whilst the one on the right is stable.

Image 1 ISTC50 test Unstable vs Stable – visual comparison



Several journal articles, posters & product promotional literature have been written regarding the new technologies available for tartrate stabilization and have quoted the Stabilab™ and its DIT% & ISTC50 results as confirmation of how effective they products work, P.Bowyer et al (2009 & 2010), D.Wollan (2010). We have found that as well as this many British retail chains like Tesco's accept very few tartrate stability tests for confirming stability and the ISTC50 is just one of the few they do accept!

Summary

We believe this to be the new standard for Tartrate stability testing as it can be used for all wine styles and gives accurate determinations no matter what stabilization process or product is used to stabilize the wine. It is relatively fast taking 5-6 hours to complete a test and provides qualitative indication of the wines instability unlike the qualitative results from the -4°C for 3 or 6 days tests. The Stabilab™ is very useful for the determination of Electrodialysis requirements, mannoprotein (Mannostab™, Laffort Australia) and CMC dosage determination via the DIT% test. Also for confirming the stability of all wines including those post treatment by the aforementioned technologies via the ISTC50 test giving a degree of comfort to those winemakers using these new technologies.

References

Stabilab™ Operating Manual. Eurodia

Handbook of Enology Volume 2: The Chemistry of Wine Stabilization and Treatments 2nd Ed. P. Ribereau-Gayon. Y. Glories. A. Maujean. D. Dubourdieu. 2006 John Wiley & Sons. Ltd. England

Wine Analysis and Production. B.W. Zoeklein. K.C. Fugelsang. B.H. Gump. F.S. Fury. 1995 Chapman & Hill. Inc New York.

Aus & NZ Grapegrower & Winemaker July 2010 issue 558 P.Bowyer. R. Marsh. C. Gouty. V. Moine. T. Battaglione. CMC: a new potassium bitartrate stabilization tool p 65-68

Aus & NZ Grapegrower & Winemaker Nov 2009 issue 550 P.Bowyer. Technical update: Mannostab – for potassium bitartrate stabilization p 105-106

Australian Technical Wine Conference 2010 Adelaide. Poster. New Tools for Potassium bitartrate precipitation inhibition – an evaluation of cellulose gum and yeast mannoproteins. P.K. Bowyer. R. Marsh. C. Gouty. V. Moine. W. Naboulet.

Aus & NZ Grapegrower & Winemaker June 2010 issue 557 D. Wollan. Electrolysis – new technology for rapid, efficient, reliable tartrate stabilization p 81-85

Rodney Marsh is the managing Director of Winechek, Western Australia and in that role provides technical advice and analytical services to the wine industry of Australia. He can be contacted on (08) 97514601 or email: rodney@winechek.com visit: www.winechek.com